


Product Datasheet

LC3C antibody (orb500569)

Description	LC3C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human LC3C (1-50/147 aa)
Target	LC3C
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16 kDa
Uniprot ID	Q9BXW4
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow
cytometry of
Hela(blue
line).Primar...