


Product Datasheet

Smac antibody (orb500576)

Description	Smac antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Smac (151-239/239 aa)
Target	Smac
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q9NR28
Dilution Range	WB=1:500-2000, ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis of SH-SY5Y cell using Smac ...