



Product Datasheet

Syndecan 2 antibody (orb500590)

Description	Syndecan 2 antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Syndecan-2 (101-201/201 aa)
Target	Syndecan 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P34741
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:K562(Human)
Cell;Primary:An...



Western blot:Lane
1:Mouse
Cerebrum;Lane ...