


Product Datasheet

Beta Amyloid 1-16 antibody (orb500591)

Description	Beta Amyloid 1-16 antibody
Species/Host	Mouse
Reactivity	Human, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide of human beta-Amyloid (1-16/42 aa)
Target	beta Amyloid 1-16
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	4.3/83 kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human U87MG Cell; Lan...