



Product Datasheet

SCLT1 antibody (orb500609)


Description	SCLT1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human SCLT1 (551-650/688 aa)
Target	SCLT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	Q96NL6
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P image of Human ovarian cancer tissue showing brown staining.



Immunohistochemical staining of paraffin-embedded tissue showing brown staining.



Immunohistochemical staining of Human tissue showing brown staining.