



Product Datasheet

CK16 antibody (orb500617)


Description	CK16 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Recombinant human CK16 protein (full length)
Target	CK16
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	51 kDa
Uniprot ID	P08779
Dilution Range	IHC-P=1:100-500, IF=1:200-800, IHC-F=1:100-500, ICC=1:100
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat uterus
using CK16 ...



IHC-P
analysis of
human skin
cancer usin...



IHC-P
analysis of
mouse skin
using CK16 ...