


Product Datasheet

HSP90 Beta antibody (orb500620)

Description	HSP90 Beta antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human HSP90 beta
Target	HSP90 beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	84 kDa
Uniprot ID	P08238
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-HSP90 b...