


Product Datasheet

Human IL-17/IL-17A ELISA Kit (orb50063)

Description	Human IL-17/IL-17A ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: I20-A155
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17504 MW
Uniprot ID	Q16552
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1
Expiration Date	Please enquire.



Human IL-17
PicoKine
ELISA Kit
standard ...