



Product Datasheet

Vimentin antibody (orb500634)


Description	Vimentin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Vimentin (151-250/466 aa)
Target	Vimentin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
CAS Number	Vimentin
Uniprot ID	P08670
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of
Hela.Primary
Antibody ...



Western
blot:U87MG
(Human)
Cell;U251(Hum...



Western
blot:U87MG(human)
Cell;Primary:A...