



Product Datasheet

TLR4 antibody (orb500648)

Description	TLR4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human TLR4 (121-200/839 aa)
Target	TLR4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90 kDa
Uniprot ID	Q9QUK6
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of K562.Primary Antibody ...



Western blot:Lane 1:H9C2 (Human) Cell;La...