


Product Datasheet

Rat Laminin ELISA Kit (orb50065)

Description	Rat Laminin ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from rat sarcoma basement membrane
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	196474 MW
Uniprot ID	F1MAN8 , D3ZN05 , A0A0G2K023 , F1LPI5 , F1MAA7 , D4A409 , F1LTF8 , D3ZQN7 , P15800 , F1LRH4 , F1M614
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat Laminin
PicoKine
ELISA Kit
standard ...