

Product Datasheet

Estrogen Related Receptor Beta antibody (orb500654)

Description	Estrogen Related Receptor Beta antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Estrogen Related Receptor beta (1-100/508 aa)
Target	Estrogen Related Receptor beta
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	095718
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.