


Product Datasheet

B7-1 antibody (orb500696)

Description	B7-1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse B7-1/CD80 (1-100/306 aa)
Target	B7-1/CD80
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	P33681
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot:Thymus
(Mouse);Primary:Anti...