

Product Datasheet

Histone H3(di Methyl K79) antibody (orb500704)

Description

Histone H3(di Methyl K79) antibody

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of di methyl K79 (DF(Di methyl K)TD)

Target

Histone H3(di methyl K79)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

MW

15 kDa

CAS Number

Histone H3.1

Uniprot ID


P68431

Dilution Range


WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:1000-5000

Expiration Date


12 months from date of receipt.



IHC-P
analysis of
rat colon
using
Histon...



IHC-P
analysis of
mouse
kidney
using His...



IHC-P
analysis of
human
colon
carcinoma
...