



Product Datasheet

GAPDH-Loading Control antibody (orb500708)


Description	GAPDH-Loading Control antibody
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombined Human GAPDH
Target	GAPDH-Loading Control
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	38 kDa
Uniprot ID	P04406
Dilution Range	WB=1:10000-1:100000, IHC-P=1:100-500, ICC=1:100
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human HEK293 cell; La...



Western blot analysis using Anti-GAPDH a...



ICC analysis of HeLa cell using GAPDH-Lo...