


Product Datasheet

Mouse M-CSF ELISA Kit (orb50071)

Description	Mouse M-CSF ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K33-E262
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA).
MW	60649 MW
Uniprot ID	P07141
Sensitivity	1 pg/ml
Dilution Range	For serum, dilution ratio is 1:4 and concentration is 510-1678 pg/ml. For plasma, dilution ratio is 1:4 and concentration is 428-1191 pg/ml.
Expiration Date	Please enquire.



Mouse M-CSF
PicoKine
ELISA Kit
standard ...