

Product Datasheet

CD34 antibody (orb500715)

Description	CD34 antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human CD34 (32-290/385 aa)
Target	CD34
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	39 kDa
Uniprot ID	P28906
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.