



Product Datasheet

Ki-67 antibody (orb500716)

Description	Ki-67 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Ki-67 (2901-3000/3256 aa)
Target	Ki-67
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	358 kDa
Uniprot ID	P46013
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=1ug/test, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of HL-60.Primary Antibody...



ICC analysis of HeLa cell using Ki-67 po...