


Product Datasheet

Human CCL22/MDC ELISA Kit (orb50072)

Description	Human CCL22/MDC ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G25-Q93
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10625 MW
Uniprot ID	000626
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml.
Expiration Date	Please enquire.



Human
CCL22/MDC
PicoKine
ELISA Kit
stand...