



Product Datasheet

TERT antibody (orb500720)


Description	TERT antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human TERT (521-620/1132 aa)
Target	TERT
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	O14746
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of Hela.Primary Antibody ...



Flow cytometry of Mouse thymus cells(blue...)



ICC analysis of A549 cell using TERT pol...