


## Product Datasheet

### StAR antibody (orb500744)

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | StAR antibody                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                               |
| <b>Reactivity</b>          | Mouse, Rat                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                         |
| <b>Tested Applications</b> | FC, IHC-P                                                                                            |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human StAR (201-300/285 aa)                            |
| <b>Target</b>              | StAR                                                                                                 |
| <b>Preservatives</b>       | Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4                        |
| <b>Form/Appearance</b>     | Liquid                                                                                               |
| <b>Concentration</b>       | 1mg/ml                                                                                               |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                      |
| <b>Note</b>                | For research use only                                                                                |
| <b>Isotype</b>             | IgG                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                           |
| <b>MW</b>                  | 32 kDa                                                                                               |
| <b>Uniprot ID</b>          | <b>P49675</b>                                                                                        |
| <b>Dilution Range</b>      | IHC-P=1:100-500, Flow-Cyt=1µg/Test, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                      |



IHC-P  
analysis of  
rat brain  
tissue using...