



Product Datasheet

HSP90 Alpha antibody (orb500766)

Description	HSP90 Alpha antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human HSP90 alpha Protein (549-692/732 aa)
Target	HSP90 alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	81 kDa
Uniprot ID	P07900
Dilution Range	WB=1:500-1000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Mouse Testis; Lane
2: ...



Western blot: HeLa
Cell (Human); NIH/3T3
Ce...