


Product Datasheet

Mouse MMP-2 ELISA Kit (orb50078)

Description	Mouse MMP-2 ELISA Kit
Reactivity	Mouse
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A30-C662
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	74102 MW
Uniprot ID	P33434
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is around 246 ng/ml.
Expiration Date	Please enquire.



Mouse MMP-2 PicoKine ELISA Kit standard ...