

Product Datasheet

Phosphotyrosine antibody (orb500799)

Description	Phosphotyrosine antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Phosphotyrosine
Target	Phosphotyrosine
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
CAS Number	21820-51-9
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.