



Product Datasheet

Vimentin antibody (orb500800)

Description	Vimentin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Vimentin (101-200/466 aa)
Target	Vimentin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
CAS Number	Vimentin
Uniprot ID	P08670
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=2µg/Test, IHC-F=1:100-500, ICC=1:100, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P analysis of human brain glioma tissue.



Western blot: HeLa (Human) Cell; Primary: Anti-Vimentin.