



Product Datasheet

Gremlin antibody (orb500806)

Description	Gremlin antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human gremlin (101-184/184 aa)
Target	Gremlin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	O60565
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Mouse
embryo
tissue us...



IHC-P
analysis of
rat rectum
using
Greml...