



Product Datasheet

Tubulin beta 3 (Neuronal Marker) antibody (orb500814)


Description	Tubulin beta 3 (Neuronal Marker) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human beta III Tubulin (401-450/450 aa)
Target	TUBB3 (Neuronal Marker)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50-55 kDa
Uniprot ID	Q13509
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100, IF=1:200-800, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of U-87MG Cells(blue).Pri...



ICC analysis of BV-2 cell using TUBB3 po...



ICC analysis of SH-SY5Y cell using TUBB3...