



Product Datasheet

Cytokeratin 8 antibody (orb500822)


Description	Cytokeratin 8 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Recombinant human CK8
Target	Cytokeratin 8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	53 kDa
Uniprot ID	P05787
Dilution Range	IHC-P=1:100-500, IF=1:200-800, Flow-Cyt=1ug/Test, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometry of MCF-7. Primary Antibody...



ICC analysis of HepG2 cell using Cytoker...



IHC-P analysis of Human colon carcinoma ...