

Product Datasheet


GPR30 antibody (orb500824)

Description	GPR30 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GPR30 (251-350/375 aa)
Target	GPR30
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	Q99527
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1µg/Test, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.




IHC-P analysis
of rat brain
using GPR30

...



IHC-P analysis
of mouse
brain using
GPR30...



Western
blot: Lane
1: Cerebrum
(Mouse); Lan...