



## Product Datasheet

### GAPDH antibody (orb500826)


|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | GAPDH antibody                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                          |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | ICC, IHC-Fr, IHC-P, WB                                                          |
| <b>Immunogen</b>           | Recombinant human GAPDH full length protein                                     |
| <b>Target</b>              | GAPDH (Loading Control)                                                         |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.           |
| <b>Form/Appearance</b>     | Liquid                                                                          |
| <b>Concentration</b>       | 1mg/ml                                                                          |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | IgG                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                      |
| <b>MW</b>                  | 38 kDa                                                                          |
| <b>Uniprot ID</b>          | <b>P04406</b>                                                                   |
| <b>Dilution Range</b>      | WB=1:10000-200000, IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |



ICC analysis of A549 cell using GAPDH (L...



IHC-P analysis of rat colon using GAPDH ...



IHC-P analysis of Human kidney using GAP...