



Product Datasheet

TLR4 antibody (orb500827)

Description	TLR4 antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH conjugated synthetic peptide derived from human TLR4 (701-800/839 aa)
Target	TLR4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90 kDa
Uniprot ID	Q9QUK6
Dilution Range	WB=1:1000-2000, IF=1:200-1000, IHC-P=1:100-500, IHC-F=1:500-1000, ELISA=1:1000-5000
Expiration Date	12 months from date of receipt.



Western
blot:Lane
1:Placenta
(Mouse);Lan...



Western blot: Lane 1: Jurkat (Human) Cell; ...