


Product Datasheet

Dopamine D2 receptor antibody (orb500835)

Description	Dopamine D2 receptor antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human DRD2 (301-400/443 aa)
Target	Dopamine D2 Receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P14416
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Cerebrum
(Mouse); Lan...