



Product Datasheet

Leptin Receptor antibody (orb500837)

Description	Leptin Receptor antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Leptin receptor (251-350/1165 aa)
Target	Leptin receptor
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	103/132 kDa
Uniprot ID	P48357
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot: Lane
1:SH-SY5Y (Human)
Cell...



Western blot:Lung
(Mouse);Primary:Anti-
L...