



Product Datasheet

ACTA2 antibody (orb500841)


Description	ACTA2 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Actin alpha
Target	alpha smooth muscle Actin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P62736
Dilution Range	WB=1:500-5000, IHC-P=1:200-800, Flow-Cyt=1ug/Test, IHC-F=1:500-1000, IF=1:500-1000
Expiration Date	12 months from date of receipt.



Flow cytometry of NIH/3T3.Primary Antibody...



IHC-P analysis of Rat stomach using alph...



IHC-P analysis of Human smooth muscle us...