


Product Datasheet

iNOS antibody (orb500853)

Description	iNOS antibody
Species/Host	Rabbit
Reactivity	Gallus, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse NOS-2 (14-100/1144 aa)
Target	iNOS
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130 kDa
Uniprot ID	P29477
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Rat liver
using iNOS
P...