



Product Datasheet

ERK1/2 antibody (orb500888)


Description	ERK1/2 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ERK1/2 (151-250/380 aa)
Target	ERK1/2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	43 kDa
Uniprot ID	Q63844
Dilution Range	WB=1:500-1000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse brain
using ERK1...



IHC-P
analysis of
mouse brain
using ERK1...



Western
blot: Lane
1:A431
(Human)
Cell;La...