


Product Datasheet

Mouse Neurotrophin-3 ELISA Kit (orb50089)

Description	Mouse Neurotrophin-3 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y139-T257
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	29588 MW
Uniprot ID	P20181
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Mouse
Neurotrophin-
3 PicoKine
ELISA Kit ...