


Product Datasheet

p63 antibody (orb500891)

Description	p63 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human P63
Target	P63
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	77 kDa
Uniprot ID	Q9H3D4
Dilution Range	Flow-Cyt=2ug/Test, WB=1:1000-5000, IHC-P=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Flow
cytometry of
A549(black
line).Prima...