



Product Datasheet

TRX antibody (orb500892)


Description	TRX antibody
Species/Host	Mouse
Reactivity	All, Human, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Thioredoxin
Target	Trx-Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P10599
Dilution Range	WB=1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Trx-Tagged fusion pr...



Western blot: Lane 1: Recombined Trx/His-t...



Western blot: Lane 1: Trx-tagged fusion pr...