



Product Datasheet

CD45 antibody (orb500902)


Description	CD45 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant human CD45
Target	CD45
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	143 kDa
Uniprot ID	P08575
Dilution Range	WB=1:500-1000, IHC-P=1:100-500, Flow-Cyt=1:50, IHC-F=1:100-500, ELISA=1:1000-5000
Expiration Date	12 months from date of receipt.



Flow cytometry of U937.Primary Antibody ...



IHC-P analysis of human tonsil using CD4...



Western blot:Lane 1:Molt-4 (Human) Cell;...