


Product Datasheet

Human OPG (TNFRSF11B) / Osteoprotegerin ELISA Kit (orb50091)

Description	Human OPG (TNFRSF11B) / Osteoprotegerin ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E22-L401
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
MW	46026 MW
Uniprot ID	O00300
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 50 pg/ml.
Expiration Date	Please enquire.



Human
OPG
PicoKine
ELISA Kit
standard
cu...