


Product Datasheet

Nestin antibody (orb500926)

Description	Nestin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Nestin (21-120/1621 aa)
Target	Nestin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	177 kDa
Uniprot ID	P48681
Dilution Range	IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1µg/Test, ELISA=1:5000-10000, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



ICC analysis of U251 cell using Nestin p...