



Product Datasheet

Caspase 8 antibody (orb500937)


Description	Caspase 8 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Caspase 8 Protein
Target	Caspase 8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	43/57 kDa
Uniprot ID	O89110
Dilution Range	WB=1:2000-5000, IHC-P=1:400-1000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat lymphoid
using Cas...



IHC-P
analysis of
mouse small
intestine ...



IHC-P
analysis of
human tonsil
using Cas...