


Product Datasheet

Leptin Receptor antibody (orb500949)

Description	Leptin Receptor antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Leptin receptor (561-660/1165 aa)
Target	Leptin receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	103/132 kDa
Uniprot ID	P48357
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1µg/Test, IHC-F=1:100-500, ICC=1:100, IF=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometry of Raji (Black).Primary A...