


Product Datasheet

Human PDGF-AB ELISA Kit (orb50095)

Description	Human PDGF-AB ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: (S87-T211) + (S82-T190)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27283 MW
Uniprot ID	P01127 , P04085
Sensitivity	2 pg/ml
Dilution Range	For serum, dilution ratio is 1:20 and concentration is 5432-34260 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 34-732 pg/ml.
Expiration Date	Please enquire.



Human
PDGF-AB
PicoKine
ELISA Kit
standar...