


Product Datasheet

Rat PDGF-AB ELISA Kit (orb50096)

Description	Rat PDGF-AB ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: (R86-K204) + (S74-T182)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27422 MW
Uniprot ID	Q05028 , P28576
Sensitivity	3 pg/ml
Dilution Range	For serum, dilution ratio is 1:5 and concentration is around 1.5 ng/ml. For plasma, dilution ratio is 1:1 and concentration is 0-46 pg/ml.
Expiration Date	Please enquire.



Rat PDGF-AB
PicoKine
ELISA Kit
standard ...