



Product Datasheet

Tau antibody (orb500964)

Description	Tau antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Tau (631-730/758 aa)
Target	Tau
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	52/79 kDa
Uniprot ID	P10636
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Lane
1:Cerebrum
(Mouse);Lan...



Western blot:U251
Cell(Human);Primary:An...