



Product Datasheet

BAX antibody (orb500966)

Description	BAX antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Bax (31-130/192 aa)
Target	Bax
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	Q07812
Dilution Range	WB=1:1000-1500, IHC-P=1:100-500, ICC=1:100, IF=1:200-500, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of HeLa cell using Bax mono...



Western blot:HeLa(Human) Cell;MCF-7(Huma...