


Product Datasheet

Human PLGF/PGF ELISA Kit (orb50097)

Description	Human PLGF/PGF ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A21-R149
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	24789 MW
Uniprot ID	P49763
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-31 pg/ml.
Expiration Date	Please enquire.



Human PLGF
PicoKine
ELISA Kit
standard c...