



Product Datasheet

Leptin Receptor antibody (orb500979)

Description	Leptin Receptor antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Leptin receptor (501-600/1165 aa)
Target	Leptin receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130 kDa
Uniprot ID	P48357
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Cerebrum (Mouse);Lung (Mous...



Western blot:Liver (Rat);Primary:Anti-Le...